

Product Features


NameTwinflax P-Preg PA11 500
Weaving.....Fabric 1/1
MatrixPolyamide 11 (100% biosourced)
Mass density.....500 g / m²
Percentage of resin.....60% by weight

Product properties

Standard followedISO 527-4 (Traction)
Test pieces**Sheets making**..... Thermocompression
.....**Preparation**Steaming 60°C for 8 hours

Tests	Value	Standard deviation
Tensile strength (MPa)	53	1
Young's traction modulus (MPa)	6540	50

Transformation parameters for thermocompression



The technical information contained in this document is considered to be accurate as of the release date. If the conditions and methods of use of the product and the information to which reference is made are beyond our control, DEHONDT COMPOSITES expressly disclaims any liability for the results obtained or arising from any use of the product. The information provided relates only to the product specifically designated and may not be applicable when this product is used in combination with other materials or in any process. The operator should thoroughly test any application before commercializing. Please see the Safety Data Sheet for information on health and safety.